

# High Frequency Muscle Massage Gun

## 數顯高頻肌肉按摩槍

### Model: MS-MC

#### Features:

- Full touch screen operation, LED display / operation simple and clear, massage, fat reduction and relaxation physical
- One key start, multi gear adjustment, long-term endurance, noise reduction design
- Deep massage, strong power, crushing pain
- Professional massage instrument, reduce muscle stiffness, and quickly restore vitality
- It is suitable for massage relaxation, recovery after exercise, warm up before training, relieve pain and meet the needs of muscle relaxation
- 全觸屏操作, LED顯示/操作簡潔明瞭, 按摩 | 減脂 | 放鬆 | 理

- 原一鍵啟動、多檔調節、持久續航、降噪設計 深層按摩、動力強勁、擊潰酸痛、 專業級按摩儀器、減輕肌肉僵硬、快速回復活力 適用於按摩放鬆、運動後恢復、訓練前熱身、緩解酸痛、滿足 多種肌肉放鬆需求



#### **Specifications**

Product Model 產品型號	MS-MC
Input voltage 輸入電壓	100-240V
Speed 轉速	3200 rpm
Color 顏色	black 黑色
Rated power 額定功率	24W
Working frequency 工作頻率	1st gear 25hz-20th gear 46hz50±10μs 1檔25Hz-20檔46Hz
Packing size 包裝尺寸	12.5*24.3*34.7cm
Package weight 包裝重量	1513g

#### Instructions

將按摩頭插入主機,底部電源開關 "ON/OFF" 打開,觸屏區輕觸開 關鍵開始工作,螢幕觸摸"+"提升檔位,"-"減小檔位 注意事項

不要把產品浸入水中,遠離液體或熱源,保持通風口沒有灰塵和碎 屑;不要拆卸設備

連續使用超過一小時,請關閉產品30分鐘,然後再次使用

充電後或使用前將充電器從電池上拔下

本產品只適合成年人,傷者請不要使用本產品

僅在身體軟組織上使用該設備,而不會引起疼痛不適。請勿將用於頭 部,身體任何堅硬或有骨頭的部位

以下人員禁用:

- >動脈瘤,出血,使用血液稀釋劑
- >帶有起搏器或除顫器的心臟病,懷孕,癌症
- >在關節內固定90天內,不得在內部固定3英寸內使用該設備。
- >敏感部位:頭部,面部,頸椎,脊柱
- >靠近骨骼且肌肉覆蓋範圍較小的位置:脛骨,腳背,手背。

www.ziglite.com

Room B3, 18/F Bonsun Industrial Building, 364-366 Sha Tsui Road, Tsuen Wan, HongKong Phone: 852-35863418 Fax: 852-35440462 info@ziglite.com