

Intelligent Moxibustion Thread Moxibustion Massager

智能灸線艾灸按摩器

Model : MS-ML

Features:

- Constant temperature control Moxa molecules overflow, directly penetrate the skin
 - Three-level constant temperature to avoid burns
 - Moxibustion function (three levels of temperature) massage function (four massage methods, 16 massage strengths)
 - Moxibustion and massage can be applied to the whole body with multi-point moxibustion
-
- 恒溫控制艾分子四溢, 直接滲透皮膚
 - 三檔恒溫避免燙傷
 - 艾灸功能(三檔溫度可調節) 按摩功能(四種按摩方式,16種按摩力度)
 - 艾灸、按摩雙管齊下全身適用多穴位可灸



Specifications

Product model 產品型號	MS-ML
Product material 產品材質	ABS+ceramic heater ABS+陶瓷加熱片
Power mode 電源方式	USB
Applicable part 適用部位	Neck 頸部位
Function 功能	Moxibustion+Massage 艾灸+按摩
Use time 使用時長	about 5-6 times 15 minutes/time 約5-6次 15分/次
Product Size 產品尺寸	580*430*20mm
Product Weight 產品重量	46g
Packing Size 包裝尺寸	205*190*50mm
Package Weight 包裝重量	330g

Instructions

- ◆ 將電極片於主機母扣連接固定, 剝掉電極片上的保護膜, 貼在身體疼痛部位(注意: 找準穴位, 貼在施灸的位置, 效果更佳)
- ◆ 長按儀器上的開關鍵, 按摩儀為開啟狀態, 短按“艾灸鍵”進入艾灸模式 短按“按摩鍵, 進入按摩模式。按“+”調節按摩力度

使用時注意事項

- ◆ 心臟病患者、劇痛性疾病患者、孕婦、血壓不正常者禁止使用。
- ◆ 切勿與其它醫療電子儀器同時使用, 如心臟起搏器、人工心肺等, 否則會導致危險。
- ◆ 請勿給兒童或無表達意識能力的人玩耍或使用, 防止燙傷。